

Contents

1	Results and discussion: autoinducer-antibiotic conjugates	2
1.1	Biological testing	2
1.1.1	Autoinducer-antibiotic conjugates	2
1.1.2	Cleavable HSL-ciprofloxacin conjugates	6
1.2	Biological testing	8
1.2.1	Antibiotic susceptibility	8
1.2.2	Quantification of biofilms	8
1.2.3	Biofilm inhibition	9
1.2.4	Biofilm dispersal	9
2	References	10

1 Results and discussion: autoinducer-antibiotic conjugates

1.1 Biological testing

1.1.1 Autoinducer-antibiotic conjugates

The eight triazoles made in ?? (see Figure 1) were tested for antibacterial and anti-biofilm activity in *P. aeruginosa* PAO1¹ and YM64.²




Figure 1

In YM64 at 5 h the HSL-ciprofloxacin conjugates **72**, **77** and **78** showed slight activity at the highest concentration, but not as much as ciprofloxacin **24**. This activity was not visible by 24 h (see Figure 3) and the compounds had no effect on biofilm formation (see Figure 4).




Figure 2: YM64 5 h.




Figure 3: YM64 24 h.




Figure 4: YM64 biofilms 24 h.

In PAO1 **78** showed similar activity to ciprofloxacin **24** at the highest concentration (see Figure 5), but not at lower concentrations. All other compounds did not show activity, and again there was no activity at 24 h or against biofilms.




Figure 5: PAO1 5 h.




Figure 6: PAO1 24 h.




Figure 7: PAO1 biofilms 24 h.

1.1.2 Cleavable HSL-ciprofloxacin conjugates

The eight cleavable HSL-ciprofloxacin conjugates, two controls and two alkynes described in ?? (see Figure 8) were tested for antibacterial and anti-biofilm activity in *P. aeruginosa* YM64.




Figure 8

Here there was more success, although the activity was still not as high as for ciprofloxacin **24**. The HSL-ciprofloxacin conjugates with *N*-(acetoxymethoxycarbonyl) linkers ($R = H$) showed activity at high concentrations. A longer linker seems to give higher activity; **93** and **94** showed activity comparable with ciprofloxacin **24** at high concentrations. Unfortunately the control **98** and alkyne **90** with *N*-(acetoxymethoxycarbonyl) linkers ($R = H$) showed higher activity than the conjugates, indicating that the HSL head wasn't contributing to the activity of the conjugates.

The conjugates with an *N*-(acetoxyethoxycarbonyl) linker ($R = Me$) did not show any activity. This suggests that they either didn't enter cells or weren't suitable substrates for esterases. The *N*-(acetoxyethoxycarbonyl) linked alkyne ($R = Me$) did show some activity, indicating that maybe it could penetrate cells more easily than the conjugates due to its lower molecular weight and/or lower polarity.




Figure 9: YM64 5 h cleavable conjugates.

1.2 Biological testing

Compounds were tested against *P. aeruginosa* PAO1[?] and YM64.[?]

The compounds were tested in PAO1, and YM64 which is missing 4 main efflux pumps. Antibacterial activity was measured by measurement of turbidity at 595 nm at 5 h and 24 h. The compounds were tested at 6 concentrations between 2 and 0.0625 μ M in LB in 96-well plates.

1.2.1 Antibiotic susceptibility

Make overnight cultures (in 10 ml LB @ 37 oC). Dilute 1/100 to make 20 ml (200 μ l culture, 19.8 ml LB) per 2 plates. Add 99 μ l diluted culture to all wells except blanks and 'Just LB' wells (96-well plate, flat-bottomed wells). Remember to make 2 if doing biofilm inhibition testing after! Add 100 μ l LB to the 'Just LB' wells. Add 1 μ l compound solutions from the master plates. Stick on adhesive cover (remove edges, DO NOT USE A LID). Record OD Place the plate in shaker @ 37 oC, 100 rpm, Record the OD every 1h for 8h, then at 16h, 24h and 48h.

1.2.2 Quantification of biofilms

Prepare two plates according to steps 1-6 of 'MIC testing'. Incubate one for 24 h and one for 48 h (this can be the one from the MIC test), both @ 37 oC, 100 rpm. Record OD Aspirate out culture (don't touch sides). Add 120 μ l water then aspirate again. Repeat x3. Add 120 μ l 0.1 Aspirate out crystal violet. Add 120 μ l water then aspirate again. Repeat x3. Add 30 Vortex plate. Read plate.

1.2.3 Biofilm inhibition

1.2.4 Biofilm dispersal

Make overnight cultures (toothpick in 10 ml LB @ 37 oC). Dilute 1/100 to make 20 ml (200 ul culture, 19.8 ml LB) per 2 plates. Add 99 ul diluted culture to all wells except the 'Just LB' wells (96-well plate, flat-bottomed wells). Make 2 plates per master plate per strain. Add 100 ul LB to the 'Just LB' wells. Shake for 24 h or 48 h @ 37 oC, 100 rpm. Record OD Add 1 ul compound solutions from the master plates. Record OD Shake for 24 h @ 37 oC, 100 rpm. Record OD Do steps 2-9 of 'Biofilm inhibition testing'.

2 References

- [1] C. K. Stover, X. Q. Pham, A. L. Erwin, S. D. Mizoguchi, P. Warrener, M. J. Hickey, F. S. L. Brinkman, W. O. Hufnagle, D. J. Kowalik, M. Lagrou, R. L. Garber, L. Goltry, E. Tolentino, Y. Yuan, L. L. Brody, S. N. Coulter, K. R. Folger, A. Kas, K. Larbig, R. Lim, K. Smith, D. Spencer, G. K. Wong, Z. Wu, I. T. Paulsen, J. Reizer, M. H. Saier, R. E. W. Hancock, S. Lory and M. V. Olson. Complete genome sequence of *Pseudomonas aeruginosa* PAO1, an opportunistic pathogen. *Nature*, 406:959–964, 2000.
- [2] Y. Morita, Y. Komori, T. Mima, T. Kuroda, T. Mizushima and T. Tsuchiya. Construction of a series of mutants lacking all of the four major *mex* operons for multidrug efflux pumps or possessing each one of the operons from *Pseudomonas aeruginosa* PAO1: MexCD-OprJ is an inducible pump. *FEMS Microbiology Letters*, 202:139–143, 2001.

Todo list